
DEVILBISS®

DeVilbiss® AutoAdjust™ LT Nasal CPAP System Instruction Guide Model 8054

CAUTION– Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
Made in U.S.A.

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 **SUNRISE
MEDICAL®**

TABLE OF CONTENTS

IEC Symbols	4
Important Safeguards	4
Travel	5
Introduction	5
What is Obstructive Sleep Apnea (OSA)	5
How Does the DeVilbiss AutoAdjust LT Nasal CPAP System Work	5
Important Parts	6
Accessories/Replacement Items	6
Set-Up	7
Operation	7
Supplemental Oxygen	10
Cleaning	
Cabinet	11
Air-Inlet Filter	11
Mask, Headgear, and Tubing	11
Maintenance	12
Troubleshooting	12
Warranty	14
Specifications	15
Declaration of Conformity	15

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IEC SYMBOLS

-  Attention, consult instruction guide.
-  Off (power: disconnection from the mains)
-  On (power: connection to the mains)
-  Type BF applied part

IMPORTANT SAFEGUARDS

When using electrical products, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

DANGER– Urgent safety information for hazards that will cause serious injury or death.

WARNING– Important safety information for hazards that might cause serious injury.

CAUTION– Information for preventing damage to the product.

NOTE– Information to which you should pay special attention.

READ ALL INSTRUCTIONS BEFORE USING.

SAVE THESE INSTRUCTIONS.

WARNING

This CPAP unit should be used only with masks recommended by the unit's manufacturer, your physician, or respiratory therapist. To avoid rebreathing of exhaled air 1) a mask should not be used unless the CPAP unit is turned on and operating properly and 2) the vent hole(s) associated with the mask should never be blocked. When the CPAP unit is turned on and functioning properly, new air from the unit flushes the exhaled air out through the mask vent hole(s). However, when the CPAP unit is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation. This warning applies to most CPAP units.

TRAVEL

This unit is equipped with a universal power supply capable of accepting line voltages of 100-240V~, 50/60Hz. No adjustments are required to the AC power connector or the line fuses to accept this entire range of voltages. If traveling to another country where the line voltage is different than the domestic supply voltage, simply obtain a line cord configured for that particular country with the appropriate international plug and substitute it for the standard plug. Operate unit as normal.

INTRODUCTION

What is Obstructive Sleep Apnea (OSA)?

Obstructive Sleep Apnea (OSA) is a condition where the upper airway relaxes during sleep causing the airway to obstruct. Apnea is the Latin word meaning "without breath." People diagnosed with sleep apnea stop breathing repeatedly during sleep, resulting in extreme exhaustion. High blood pressure, cardiac enlargement, and cardiac arrhythmias are common, as are shortness of breath and loud snoring. In a 1993 report, the National Commission on Sleep Disorders Research stated that as many as 18 million Americans may have signs and symptoms of sleep apnea.

How does the DeVilbiss AutoAdjust LT nasal CPAP system work?

The DeVilbiss AutoAdjust LT nasal CPAP instrument provides non-invasive, continuous air pressure to your upper airway, acting as an air "splint" to keep the airway passages open during sleep. This reduces and eliminates the obstructive events, allowing uninterrupted sleep. The DeVilbiss AutoAdjust LT is designed to provide this therapy to you quietly and comfortably so that your sleep is restorative and refreshing. The DeVilbiss AutoAdjust LT nasal CPAP system offers automatic adjustment of air pressure by sensing your air flow and compensating to meet your needs. AutoAdjust senses snoring and changes in breathing flow and continually adjusts CPAP pressures in relation to these changes. AutoAdjust increases CPAP pressures when snoring or specific breath flow patterns are detected, and decreases the CPAP pressure during periods of sleep when snoring or abnormal breath flows are not detected. The pressure profile changes are gradual to prevent arousals of sleep architecture, but are timely to respond to respiratory events requiring pressure adjustments.

NOTE– This product has been designed for treatment of obstructive sleep apnea, specifically for use in the adult patient population.

IMPORTANT PARTS

Figures 1 & 2

1. A/D (Activate/Deactivate) Button
2. Outlet Port
3. Carrying Case
4. Air-Inlet Filter
5. Power Switch
6. AC Power Connector
7. DIN connector
8. Power Cord (assure proper configuration for outlet)
9. Service Required LED


Figure 1

ACCESSORIES/REPLACEMENT ITEMS

Figure 3

- | | |
|--|-----------|
| 1. CPAP Humidifier | 900HC105 |
| Distributed by Sunrise Medical in the U.S. only. | |
| 2. Nasal Masks | |
| Small..... | 9351D-668 |
| Medium Small..... | 9351D-667 |
| Medium | 9351D-669 |
| Large..... | 9351D-670 |
| 3. Air-Inlet Filter (4/pk) | 8000D-602 |
| Fine Particle Air-Inlet Filter | |
| (not shown) (4/pk) | |
| | 8000D-603 |
| 4. Oxygen Adapter | 7353D-601 |
| 5. DeVilbiss Seal-Ring® | |
| Small..... | 9351D-658 |
| Medium | 9351D-659 |
| Large..... | 9351D-660 |
| 6. 6' 22 mm Tubing | 7351D-616 |
| 7. Strap-Style Headgear | 7351D-617 |
| 8. Cap-Style Headgear | 7351D-627 |
| 9. DeVilbiss Humidifier..... | 9100D |
| 10. Serenity™ Mask | 9352D |
| Replacement Cushion..... | 9352D-601 |
| 15mm Tubing Swivel..... | 9352D-602 |
| Mask Adapt O ₂ Port | 9352D-603 |
| Replacement Pads and Vent | 9352D-604 |
| Mask Headgear Strap/Mesh | 9352D-619 |
| Serenity Mask w/out Headgear | 9352D-621 |


Figure 2


Figure 3

Ask your Sunrise Medical provider for information about other DeVilbiss mask styles and sizes.

SET-UP

WARNING

This is not a life support device and is intended to be used on spontaneously breathing adults.

NOTE– Refer to *Important Parts* for part identification.

- Place unit on level, sturdy surface away from heating vents, radiators, drapes, and bedspreads so that air intake is not restricted.
- Insert the appropriate power cord into the AC power connector (**Figure 2**).
- Make sure the power switch is in the “Off” position. Plug the power cord into an appropriate wall outlet.
- Connect the 6-foot length of 22 mm diameter tubing to the outlet port on the front of the unit **Figure 4**. This is the port from which CPAP air supply is provided.
- Connect the opposite end of the 22 mm tubing to the mask and headgear assembly.

NOTE– Excess facial oils will deteriorate mask material. To prolong mask life, wash face nightly before wearing mask.

NOTE– Some tubing is prone to kinking. Flow can be blocked causing apnea episodes to resume. Use only 22 mm diameter smooth bore-type, kink-resistant tubing.

- Position mask loosely on face. Adjust headgear straps so mask is snug and comfortable **Figure 5**.

OPERATION

Step 1

- Press power switch (rear of unit) to the “On” position. The blower will operate.

Step 2

- Your DeVilbiss AutoAdjust LT nasal CPAP device will automatically adjust the CPAP pressure up and down throughout the night based on your needs. If you wake up and the pressure seems too high, simply press the A/D button on the front of the unit and the pressure will return to a low pressure operating mode **Figure 6**.

Step 3

- Breathe into the mask and go to sleep.


Figure 4

Figure 5



Figure 6

OPERATION

The DeVilbiss AutoAdjust LT nasal CPAP system, model #8054, offers automatic adjustment of the administered CPAP pressure. The device is equipped with a pneumotachometer which monitors the patient's air flow signal for purposes of identifying any flow restrictions and snoring episodes often associated with respiratory events. The DeVilbiss AutoAdjust LT nasal CPAP device has been designed to minimize these flow restrictions which often lead to patient arousal and fragmented sleep by automatically increasing the CPAP pressure in the presence of respiratory events. A proprietary pressure adjustment control algorithm has been developed by Sunrise Medical which is based around event density and event type.

The pressure profile changes are gradual so as to prevent arousals of sleep architecture, but are timely to respond to respiratory events requiring pressure adjustments. The device also reduces airway pressure in an attempt to minimize the average supplied pressure administered throughout the night. Refer to **Diagram 1** of a typical pressure profile administered to a patient and associated operating states.

- DIAGRAM 1
- A) Power switch turned on.
 - A-B) AutoAdjust is in low pressure pause state, providing low level of flow and CPAP pressure.
 - B) Patient breathing is detected, pressure is increased to the lower pressure limit prescription setting.
 - C-D) Patient falls asleep.
 - D-E) Patient begins to have respiratory events. Unit senses the events and adjusts pressure.
 - E-F) Pressure returns to the lower limit as events become less frequent. (**Note:** Pressure cannot be less than the lower pressure limit unless the A/D button is pressed to reach the low pressure pause mode.)
 - F-G) Respiratory events resume and the CPAP pressure is automatically adjusted to reduce the number of events.
 - G) Patient presses the A/D button on the front of the unit, causing the pressure to drop to low pressure pause level. The low pressure pause mode of operation provides a low level of flow and pressure (approximately 3 cm H₂O). If the patient wakes up with a high CPAP pressure, this allows the patient to fall back asleep more comfortably.
 - G-H) Patient falls back asleep and begins to have respiratory events. The AutoAdjust unit detects the events and automatically adjusts the CPAP pressure level.
 - H) CPAP pressure has been automatically increased and has reached the upper pressure limit prescription setting for the product. (**Note:** CPAP pressure cannot go over the upper limit regardless of patient events.)
 - I) A large leak in the system, usually at the patient mask, has been detected by the AutoAdjust CPAP and the unit has reduced the pressure to the low pressure pause level.
 - J-K) Breathing signals have been sensed at the mask and normal operation of the unit has resumed. When breathing is detected, the pressure is raised to 80% of the pressure level that the unit was operating at when the leak was sensed.
 - K) The patient has awakened and has pressed the A/D button on the front of the unit to reduce the pressure to the low pressure pause level.
 - L) The patient switches the unit off.

SUPPLEMENTAL OXYGEN

Supplemental oxygen may be prescribed with nasal CPAP for treatment of obstructive sleep apnea.

CAUTION– Oxygen is a prescription gas and should only be administered under the supervision of a physician.

⚠ WARNING

Oxygen supports combustion. To avoid possible physical injury, do not smoke while using your CPAP unit with supplemental oxygen. Do not use near hot, sparking objects or open sources of flame.

⚠ WARNING

Always turn the CPAP unit on prior to turning on the oxygen flow and turn off the oxygen before turning the unit off.

1. Activate power switch and connect the optional oxygen adapter (part #7353D-601) to the outlet port of the unit. Always turn CPAP unit on before the oxygen source and turn the oxygen source off before turning off the CPAP unit.
2. Connect one end of the 22 mm diameter tubing to the open end of the adapter and the other end to the CPAP mask.
3. Oxygen may be injected into the air supply by connecting one end of standard oxygen tubing to the oxygen adapter and the other end to the oxygen source **Figure 7**.

⚠ WARNING

The flow setting of the oxygen source shall be specified by a physician.

⚠ WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

NOTE– Equipment is suitable for use with oxygen if prescribed by your physician.


Figure 7

CLEANING

⚠ WARNING

To avoid electric shock, always unplug power cord from wall outlet power source.

CAUTION– Never rinse or place the unit in water. Never allow liquids to get into or around any of the ports, switches, or air filter; doing so will result in product damage.

Cabinet

1. Unplug the power cord from the wall outlet.
2. Wipe the cabinet with a clean, damp cloth and mild dishwashing detergent every few days to keep dust free **Figure 8**. Do not use alcohol or any alcohol-based solutions to clean the unit. These cleaning solutions may damage the plastic enclosure. **Allow the unit to dry completely before plugging the unit into a power source.**

Air-Inlet Filter

1. The air-inlet filter should be checked daily. The filter is located in the rear of the unit **Figure 9**.
2. For weekly cleaning, remove the dirty foam filter from the frame **Figure 10**.
3. Wash the filter in a solution of warm water and dishwashing detergent and rinse with tap water. Allow the filter to dry completely. If the filter (part #8000D-602) is damaged, contact your home medical equipment provider for a replacement.
4. Insert clean filter into housing **Figure 11**.

Mask, Headgear, and Tubing

CAUTION– Do not place mask, headgear, or tubing in an automatic dishwasher. Some automatic dishwashing detergents contain chemicals that may damage these items. Also, do not clean with alcohol or any alcohol-based solutions as damage will result.

1. Remove the tubing and headgear from the mask.
2. Each of these items may be washed as often as needed by using a solution of warm water and dishwashing detergent. Wash item individually, rinse with clean, warm tap water and allow to air dry before using. Hang the tubing vertically to allow the water to drain and dry completely.


Figure 8


Figure 9


Figure 10


Figure 11

MAINTENANCE

⚠ DANGER

Electric shock hazard. Do not attempt to open or remove cabinet, there are no user-serviceable internal components. If service is required, return it to a qualified Sunrise Medical provider or authorized service center. Opening or tampering with the unit will void warranty.

TROUBLESHOOTING

NOTE– Your DeVilbiss AutoAdjust LT nasal CPAP system contains no user-serviceable parts. If you believe your unit is not working properly, BEFORE YOU RETURN IT TO THE HOME MEDICAL EQUIPMENT PROVIDER WHERE YOU PURCHASED IT OR TO SUNRISE MEDICAL, please take a few moments to check for these possible causes:

BASE UNIT		
Symptoms	Possible Causes	Remedies
Unit does not turn on.	<ol style="list-style-type: none"> Unit not connected to a live power source. Power cord not attached properly. 	<ol style="list-style-type: none"> Plug unit into a live receptacle. Firmly attach power cord.
No air flow emitted from flow generator when power is turned on, motor runs.	<ol style="list-style-type: none"> Air filter may be occluded by dirt or bed clothing. 	<ol style="list-style-type: none"> Clean or replace the air filter. Move unit away from loose material.
Symptoms of OSA have recurred.	<ol style="list-style-type: none"> Air filter may be occluded. Incorrect air pressure setting. Sleep apnea condition has changed. 	<ol style="list-style-type: none"> Clean or replace the air filter and relocate unit away from drapes, etc. Contact your Sunrise Medical provider or authorized service center. Call your physician; you may need to be retested.
Skin becomes irritated where mask touches skin.	<ol style="list-style-type: none"> Headgear is too tight. Incorrect mask size. Allergic reaction to the mask material. Worn mask. 	<ol style="list-style-type: none"> Loosen the headgear. Contact your Sunrise Medical provider. Call your physician. Replace mask.
Dryness of throat or nose.	<ol style="list-style-type: none"> Inadequate room humidity. Room temperature too high. 	<ol style="list-style-type: none"> Increase room humidity or add a CPAP humidifier. Decrease room temperature.
Air from flow generator seems warm.	<ol style="list-style-type: none"> Air filter may be occluded. Air filter and/or outlet port may be obstructed. Unit located too close to a heating vent or radiator. 	<ol style="list-style-type: none"> Clean or replace the air filter. Move the unit away from obstruction. Move the unit away from heat source.
Nasal, sinus, ear pain, or runny nose.	<ol style="list-style-type: none"> Reaction to air flow pressure, sinus infection, or nasal congestion. 	<ol style="list-style-type: none"> Discontinue use. Contact your physician.
Yellow light on back panel of unit is activated.	<ol style="list-style-type: none"> Internal component failure. 	<ol style="list-style-type: none"> Contact your Sunrise Medical provider or authorized service center.
Aroused from sleep by a sudden drop in pressure followed by a sudden rise.	<ol style="list-style-type: none"> Device measured an extremely high leak level in the mask. 	<ol style="list-style-type: none"> Readjust mask to seal leak or replace old, worn mask with a new one.

WARRANTY

DeVilbiss AutoAdjust LT Model 8054 nasal CPAP system is warranted to be free from defective workmanship and material for a period of two years from date of purchase. Any defective part(s) will be repaired or replaced at Sunrise Medical's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your Sunrise Medical provider or Sunrise Medical Service Department for instructions:

U.S.A. 800-333-4000, 814-443-4881

Canada 905-660-2459

Europe 44-138-444-6688

NOTE—This warranty does not cover providing a loaner unit, compensating for costs incurred in rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

SPECIFICATIONS

Size4.0"H x 7.7"W x 10.3"D
(10.2 cm x 19.6 cm x 26.2 cm)

Weight3.6 lbs. (1.7 kg)

Electrical Requirements.....100-240V~, 50/60 Hz

Power Consumption60 watts maximum

CPAP Pressure Range3-20 cm H2O

Operating Temperature Range41°-97°F (5° - 36°C)

Operating Humidity Range95% RH (noncondensing)

Operating Atmospheric ConditionsSea level to 10,000 feet

Storage & Transportation Temperature Range-4 to 140°F (-20 - 60°C)

Storage & Transportation Humidity Range95% RH (noncondensing)

WarrantyTwo-year limited

Equipment classification with respect to protection from electric shock:....Class II

Degree of protection from electric shock.....Type BF

Degree of protection against ingress of liquidsNone

Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

Mode of operationContinuous.

DECLARATION OF CONFORMITY

Manufacturer: Sunrise Medical

Address: Respiratory Products Division
100 DeVilbiss Drive
Somerset, Pennsylvania 15501
USA

Product Designation: CPAP

Type, Model:
DeVilbiss® AutoAdjust™ LT 8054

We herewith declare that the above-mentioned product complies with the requirements of EC Directive 93/42/EEC and the following:

**Quality System
Standards Applied:** IS09001/EN46001

Notified Body RWTÜV

MDD Annex II Applied

**Safety Standards
Applied:** IEC 601-1/EN60601-1

EMC Compliance to:

EMC Directive 89/336/EEC
IEC 801-2 through 5
CISPR 11/EN55011

Place: Somerset, PA

Signature: 

Name: Michael Chellson

Title: Director of quality and regulatory affairs

European Contact:

Sunrise Medical Ltd.
Sunrise Business Park
High Street
Wollaston, West Midlands
DY8 4PS
ENGLAND
44-138-444-6688

®Registered U.S. Patent and Trademark Office and other countries.

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©2001, Sunrise Medical 6.01
Sunrise Medical Pty. Limited