
DEVILBISS®

DeVilbiss® Nasal CPAP System Instruction Guide 9000 & 9001 Series

CAUTION—Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. Made in U.S.A.

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 **SUNRISE
MEDICAL®**

TABLE OF CONTENTS

IEC Symbols	3
Important Safeguards	4
Travel	4
Introduction	5
Important Parts	6
Accessories/Replacement Items	6
Set-Up for 9000/9001	7
Operation for 9000	7
Operation for 9001	9
DC Operation for 9001	11
Supplemental Oxygen	12
Cleaning	13
Cabinet	13
Air-Inlet Filter	13
Mask, Headgear and Tubing	14
Maintenance	14
Troubleshooting	15
Specifications	16
Warranty	20

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
IEC Symbols

-  ON / OFF- turns blower on and off
-  Select Next Item on LCD display
-  Increase value shown on LCD display
-  Decrease value shown on LCD display
-  Start Delay- starts the delay cycle, if prescribed
-  Attention, consult instruction guide
-  Type BF equipment- applied part
-  Class II electrical protection- double insulated
-  Data port
-  DC Input- 12 VDC, center pin is positive
-  100 – 240 Volts AC input 50/60 Hz

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IMPORTANT SAFEGUARDS

When using electrical products, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

▲ DANGER Urgent safety information for hazards that will cause serious injury or death.

▲ WARNING Important safety information for hazards that might cause serious injury.

CAUTION Information for preventing damage to the product.

NOTE Information to which you should pay special attention.

READ ALL INSTRUCTIONS BEFORE USING.

SAVE THESE INSTRUCTIONS.

▲ WARNING

This CPAP unit should be used only with masks recommended by the unit's manufacturer, your physician, or respiratory therapist. To avoid rebreathing of exhaled air 1) a mask should not be used unless the CPAP unit is turned on and operating properly, and 2) the vent hole(s) associated with the mask should never be blocked. When the CPAP unit is turned on and functioning properly, new air from the unit flushes the exhaled air out through the mask vent hole(s). However, when the CPAP unit is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation. This warning applies to most CPAP units.

CAUTION— At low CPAP pressures, the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.

TRAVEL

This unit is equipped with a universal power supply capable of accepting line voltages of 100-240V~, 50/60Hz (9000/9001) and 12V (9001 only). No adjustments are required to the AC power connector or the line fuses to accept this entire range of voltages. If traveling to another country where the line voltage is different than the domestic supply voltage, simply obtain a line cord configured for that particular country with the appropriate international plug and substitute it for the standard plug. Operate unit as normal.

To find out which model you have either look at the label on the bottom of the unit or look at the outlet port color **Figure 1**. Model 9000 has a grey outlet port. Model 9001 has a black outlet port.


Figure 1

▲ DANGER

Electric shock hazard - Do not use while bathing.

Electric shock hazard - Do not place or drop into water or any other liquid.

NOTE—Your home medical equipment provider has set your prescription pressure based on the altitude in your area. Model 9000 needs to be manually changed for altitude compensation, 9001 has onboard automatic altitude compensation. If an adjustment is necessary, be sure to readjust the pressure upon return.

Set the altitude range on the LCD menu based on the chart below. Model 9000 can be adjusted by pressing the select key to advance to the altitude setting. The proper altitude can then be set by pressing the up or down key.

Altitude	LCD Symbol
0-2500	1 ■
2500-5000	2 ■■
5000-7500	3 ■■■
7500-9000	4 ■■■■

INTRODUCTION

What is Obstructive Sleep Apnea (OSA)?

Obstructive Sleep Apnea (OSA) is a condition where the upper airway relaxes during sleep causing the airway to obstruct. Apnea is the Latin word meaning “without breath.” People diagnosed with sleep apnea stop breathing repeatedly during sleep, resulting in extreme exhaustion. High blood pressure, cardiac enlargement, and cardiac arrhythmias are common, as are shortness of breath and loud snoring. In a 1993 report, the National Commission on Sleep Disorders Research stated that as many as 18 million Americans may have signs and symptoms of sleep apnea.

How does the nasal CPAP system work?

The CPAP instrument provides non-invasive continuous air pressure to your upper airway, acting as an air “splint” to keep the airway passages open during sleep. This reduces and eliminates the obstructive events, allowing uninterrupted sleep. The DeVilbiss 9000/9001 Series is designed to provide this therapy to you quietly and comfortably so that your sleep is restorative and refreshing.

NOTE—This product has been designed for the treatment of obstructive sleep apnea, specifically for use in the adult patient population.

IMPORTANT PARTS - FIGURES 2 & 3

1. Keypad/LCD Display (refer to **Figure 5** for LCD symbol definition)
2. Outlet Port
3. Carrying Case
4. Air-Inlet Filter
5. DC input connector (model 9001 only)
6. AC Power Connector
7. DIN connector
8. Power Cord (assure proper configuration for outlet)


Figure 2

ACCESSORIES/REPLACEMENT ITEMS—FIGURE 4

- | | |
|--|-----------|
| 1. CPAP Humidifier | 9100D |
| F & P Humidifier (not shown) | HC105 |
| 2. DeVilbiss Serenity™ Mask | 9352D |
| (One size fits all. Includes headgear.) | |
| Nasal CPAP Masks (not shown) | |
| Small | 9351D-668 |
| Medium | 9351D-669 |
| Medium/Small | 9351D-667 |
| Large | 9351D-670 |
| 3. Air-Inlet Filter(4/pk) | 8000D-602 |
| 4. Optional Fine Particle Filter | 8000D-603 |
| 5. Oxygen Adapter | 7353D-601 |
| 6. DeVilbiss Seal-Ring™ (4/pk) (not shown) | |
| Small | 9351D-658 |
| Medium or Medium/Small | 9351D-659 |
| Large | 9351D-660 |
| 7. 6' 22 mm Tubing | 7351D-616 |
| 8. Strap-Style Headgear (not shown) | 7351D-617 |
| 9. Cap-Style Headgear (not shown) | 7351D-627 |
| 10. Optional 12V DC power cord (9001 only) | 7304D-619 |


Figure 3


Figure 4

Ask your Sunrise Medical provider for information about other DeVilbiss mask styles and sizes.

SET-UP FOR 9000/9001**▲ WARNING**

This is not a life support device and is intended to be used on spontaneously breathing adults.

NOTE—Refer to *Important Parts* for part identification.

- Place unit on level, sturdy surface away from heating vents, radiators, drapes and bedspreads so that air intake is not restricted.
 - Insert the appropriate power cord into the AC power connector (**Figure 3**).
 - Plug the power cord into an appropriate wall outlet. The display will show *Set Pres* and CPAP will increase to the prescribed pressure.
 - Refer to **Figure 5** for LCD symbol definitions.
 - Connect the 6-foot length of 22 mm diameter tubing to the outlet port on the front of the unit (**Figure 6**). This is the port from which CPAP air supply is provided.
 - Connect the opposite end of the 22 mm tubing to the mask and headgear assembly.
- NOTE**—Excess facial oils will deteriorate mask material. To prolong mask life, wash face nightly before wearing mask.
- NOTE**—Some tubing is prone to kinking. Flow can be blocked causing apnea episodes to resume. Use only 22 mm diameter smooth bore-type, kink-resistant tubing.
- Position mask loosely on face. Adjust headgear straps so mask is snug and comfortable (**Figure 7**).

OPERATION OF THE 9000**CO₂ Rebreathing**

CO₂ rebreathing is minimized during use by providing 3cmH₂O minimum pressure, which creates airflow through mask vent holes. It is important that mask vent holes are not blocked.


Figure 5


Figure 6


Figure 7

During a single fault, the blower is maintained within spec or at less than the maximum limited pressure. For a firmware detected fault, the blower is set to a fail safe mode in which the blower is turned off and the detected error code is shown on the LCD display. CO₂ rebreathing is minimized by keeping the mask vent holes open at all times.

1. Plug CPAP into AC power using appropriate power cord. Will work on 100-240 volts AC (50 or 60 Hz) with no changes except for the end of the power cord where it plugs into the wall.
2. The display will show *Set Pres* and the CPAP will increase to the prescribed CPAP pressure.
3. The display shows the CPAP pressure prescribed by the doctor.

NOTE—Your physician has chosen the appropriate delay time (0, 10, 20, 30 or 45 minutes) that has been preset by your Sunrise Medical provider or clinician.

4. Press the Delay button (**Figure 8**) to start the pressure delay cycle if prescribed. Pressure will decrease to the Delay Pressure setting and return to the prescribed Set Pressure after the Delay Time has expired.
5. When you awake in the morning, press the ON/OFF button to turn the blower OFF (**Figure 9**).

Model 9000 Altitude Selection

Model 9000 has manual altitude selection. This allows the patient to correct the CPAP pressure when traveling.

To correct the CPAP pressure for altitude changes:

1. Plug CPAP into AC power using appropriate power cord. Will work on 100 to 240 volts AC (50 or 60 Hz) with no changes except for the end of the power cord where it plugs into the wall.
2. The display will show *Set Pres* and the CPAP will increase to the prescribed CPAP pressure.
3. Press the Select button on the keypad until the LCD display shows *Altitude*.

Altitude	LCD Symbol
0-2500	1 ■
2500-5000	2 ■■
5000-7500	3 ■■■
7500-9000	4 ■■■■


Figure 8


Figure 9

4. Use the UP or DOWN keys on the keypad to change the altitude setting.
5. Operate the CPAP normally.
6. Remember to return *Altitude* to its original setting upon returning home.

Model 9000 LCD Display Window Definitions

Pressing the SELECT button changes the display to show information about the CPAP. Press the UP or DOWN button to change the displayed setting. (The patient can only change the Altitude and Language.)

- Set Pressure— CPAP pressure (cmH₂O) prescribed by your doctor
- Altitude Range— Patient selectable, corrects CPAP pressure for changes in altitude while traveling. See Operating Instructions for ranges and use.
- Delay Time— Time for the CPAP to reach the prescription pressure after pressing the Delay button.
- Delay Pressure— Pressure delivered during the delay cycle.
- Compliance Hours— Number of hours recorded while breathing was present.
- Unit Hours— Number of hours recorded while CPAP blower was ON.
- Firmware Version— Identification for firmware installed in the CPAP.
- Serial Number— Serial number of the CPAP unit, same as on the label.
- Language— Language selected for LCD display messages.
- Off- Blower is OFF AC or DC power is connected.
- Error- An error has been detected by the firmware, the CPAP is in a fail safe mode.

OPERATION OF THE 9001

CO₂ Rebreathing

CO₂ rebreathing is minimized during use by providing 3cmH₂O minimum pressure, which creates airflow through mask vent holes. It is important that mask vent holes are not blocked.

During a single fault, the blower is maintained within spec or at less than the maximum limited pressure. For a firmware detected fault, the blower is set to a fail safe mode in which the blower is turned off and the detected error code is shown on the LCD display. CO₂ rebreathing is minimized by keeping the mask vent holes open at all times.

1. Plug CPAP into AC power using appropriate power cord. Will work on 100 to 240 volts AC (50 to 60 Hz) with no changes except for the end of the power cord where it plugs in the wall. Model 9001 will also operate on 12 VDC.
2. The display will show *Set Pres* and the CPAP will increase to the prescribed CPAP pressure.
3. The display shows the CPAP pressure prescribed by the doctor.
4. Press the Delay button to start the pressure delay cycle, if prescribed. Pressure will decrease to the Delay Pressure setting and return to the prescribed Set pressure after the Delay Time has expired (**Figure 8**).
5. When you awake in the morning, press the ON/OFF button to turn the blower OFF (**Figure 9**).

Model 9001 also has automatic start up and automatic shut off features. To start the blower, connect the 6 foot hose and mask, then breathe into the mask. The blower will start after the 2nd or 3rd breath. To turn the blower OFF, remove the mask. After about one minute of mask off, the blower will turn OFF.

Model 9001 LCD Display Window Definitions

Pressing the SELECT button changes the display to show information about the CPAP. Press the UP or DOWN button to change the displayed setting. (The patient can only change the Language.)

- Set Pressure– CPAP pressure (cmH₂O) prescribed by your doctor
- Delay Time– Time for the CPAP to reach the prescription pressure after pressing the Delay button.
- Delay Pressure– Pressure delivered during the delay cycle.
- Compliance Hours– Number of hours recorded while breathing was present.
- Unit Hours– Number of hours recorded while CPAP blower was ON.
- Firmware Version– Identification for firmware installed in the CPAP.
- Serial Number– Serial number of the CPAP unit, same as on the label.
- Language– Language selected for LCD display messages.
- Off- Blower is OFF AC or DC power is connected.
- Error- An error has been detected by the firmware, the CPAP is in a fail safe mode.


DC Operation of the 9001

The DeVilbiss CPAP Model 9001 can be operated from a 12V DC power source in a motor home, recreational vehicle or boat. Connect the DeVilbiss CPAP to a 12V cigarette lighter in a vehicle with a negative ground using DeVilbiss adapter plug part #7304D-619. There are no switches or fuses to change for DC operation. Do not use the CPAP in a moving vehicle.

CAUTION– If this device is used with a 12VDC power source, loss of power will result in loss of CPAP therapy; there is no internal battery. If the 12VDC power supply is an automobile or truck battery, the battery should be professionally tested to make sure the battery is properly charged. A weak battery can result in insufficient voltage to operate the CPAP and loss of CPAP therapy. Do not operate the unit while the vehicle is running.

CAUTION– Use only DeVilbiss DC power cord part #7304D-619. Use of any other cord may cause damage to the CPAP or the vehicle electrical system.

CAUTION– The CPAP will operate for 8 hours on a car battery. The battery should be professionally tested before using the CPAP overnight to make sure that the battery is in good condition and is fully charged.

SUPPLEMENTAL OXYGEN

Supplemental oxygen may be prescribed with nasal CPAP for treatment of obstructive sleep apnea.

CAUTION—Oxygen is a prescription gas and should only be administered under the supervision of a physician.

▲ WARNING

Oxygen supports combustion. To avoid possible physical injury, do not smoke while using your CPAP unit with supplemental oxygen. Do not use near hot, sparking objects or open sources of flame.

▲ WARNING

Always turn the CPAP unit on prior to turning on the oxygen flow and turn off the oxygen before turning the unit off.

1. Activate power switch and connect the optional oxygen adapter (part #7353D-601) to the outlet port of the unit. Always turn CPAP unit on before the oxygen source and turn the oxygen source off before turning off the CPAP unit.
2. Connect one end of the 22 mm diameter tubing to the open end of the adapter and the other end to the CPAP mask.
3. Oxygen may be injected into the air supply by connecting one end of standard oxygen tubing to the oxygen adapter and the other end to the oxygen source (Figure 10).

▲ WARNING

The flow setting of the oxygen source shall be specified by a physician.

▲ WARNING

Equipment is not suitable in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

NOTE—Equipment is suitable for use with oxygen if prescribed by your physician.


Figure 10

NOTE—At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on the pressure setting, patient breathing pattern, mask selection and leak rate. This applies to most types of CPAP machines.

CLEANING**▲ WARNING**

To avoid electric shock, always unplug power cord from wall outlet power source.

CAUTION—Never rinse or place the unit in water. Never allow liquids to get into or around any of the ports, switches, or air filter; doing so will result in product damage.

Cabinet

1. Unplug the power cord from the wall outlet.
2. Wipe the cabinet with a clean, damp cloth and mild dishwashing detergent every few days to keep dust free. Do not use alcohol or any alcohol-based solutions to clean the unit. These cleaning solutions may damage the plastic enclosure. **Allow the unit to dry completely before plugging the unit into a power source.**

Air-Inlet Filter

1. The air-inlet filter should be checked daily. The filter is located in the rear of the unit (Figure 11).
2. For weekly cleaning, remove the dirty foam filter from the frame.
3. Wash the filter in a solution of warm water and dishwashing detergent and rinse with tap water. Allow the filter to dry completely. If the filter (part #8000D-602) is damaged, contact your home medical equipment provider for a replacement.
4. Insert clean filter into housing.


Figure 11

Mask, Headgear, and Tubing

CAUTION—Do not place mask, headgear or tubing in an automatic dishwasher. Some automatic dishwashing detergents contain chemicals that may damage these items. Also, do not clean with alcohol or any alcohol-based solutions as damage will result.

1. Remove the tubing and headgear from the mask.
2. Each of these items may be washed as often as needed by using a solution of warm water and dishwashing detergent. Wash item individually; rinse with clean, warm tap water and allow to air dry before using. Hang the tubing vertically to allow the water to drain and dry completely.

MAINTENANCE

▲ DANGER

Electric shock hazard. Do not attempt to open or remove cabinet; there are no user-serviceable internal components. If service is required, return it to a qualified Sunrise Medical provider or authorized service center. Opening or tampering with the unit will void warranty.

Maintenance Interval

Pressure— Check pressure once per year with a calibrated gauge, calibrate unit if:

Model 9000 error > ± 2 cmH₂O

Model 9001 error > ± 1 cmH₂O

Standard Filter— Check filter daily, clean if necessary. Replace every 6 months or if damaged.

Optional Fine Particle Filter— Check filter daily, replace if dirty. Replace every 30 days or if damaged.

TROUBLESHOOTING

NOTE—DeVilbiss 9000/9001 Series contains no user-serviceable parts. If you believe your unit is not working properly, BEFORE YOU RETURN IT TO THE HOME MEDICAL EQUIPMENT PROVIDER WHERE YOU PURCHASED IT OR TO SUNRISE MEDICAL, please take a few moments to check for these possible causes:

BASE UNIT		
Symptoms	Possible Causes	Remedies
Unit does not turn on.	<ol style="list-style-type: none"> 1. Unit not connected to a live power source. 2. Power cord not attached properly. 	<ol style="list-style-type: none"> 1. Plug unit into a live receptacle. 2. Firmly attach power cord.
No air flow emitted from flow generator when power is turned on, motor runs.	<ol style="list-style-type: none"> 1. Air filter may be occluded by dirt or bed clothing. 	<ol style="list-style-type: none"> 1a. Clean or replace the air filter. 1b. Move unit away from loose material.
Symptoms of OSA have recurred.	<ol style="list-style-type: none"> 1. Air filter may be occluded. 2. Incorrect air pressure setting. 3. Sleep apnea condition has changed. 	<ol style="list-style-type: none"> 1. Clean or replace the air filter and relocate unit away from drapes, etc. 2. Contact your Sunrise Medical provider or authorized service center. 3. Call your physician; you may need to be retested.
Skin becomes irritated where mask touches skin.	<ol style="list-style-type: none"> 1. Headgear is too tight. 2. Incorrect mask size. 3. Allergic reaction to the mask material. 4. Worn mask. 	<ol style="list-style-type: none"> 1. Loosen the headgear. 2. Contact your Sunrise Medical provider. 3. Call your physician. 4. Replace mask.
Dryness of throat or nose.	<ol style="list-style-type: none"> 1. Inadequate room humidity. 2. Room temperature too high. 	<ol style="list-style-type: none"> 1. Increase room humidity or add a CPAP humidifier; see page 6. 2. Decrease room temperature.
Air from flow generator seems warm.	<ol style="list-style-type: none"> 1. Air filter may be occluded. 2. Air filter and/or outlet port may be obstructed. 3. Unit located too close to a heat register or radiator. 	<ol style="list-style-type: none"> 1. Clean or replace the air filter. 2. Move the unit away from obstruction. 3. Move the unit away from heat source.
Nasal, sinus, ear pain, or runny nose.	<ol style="list-style-type: none"> 1. Reaction to air flow pressure, sinus infection, or nasal congestion. 	<ol style="list-style-type: none"> 1. Discontinue use. Contact your physician.

SPECIFICATIONS

Technical Specifications

Size4.0"H x 7.5"W x 10.8"D
 (10.2 cm x 19.6 cm x 26.2 cm)

Weight.....3.6 lbs. (1.7 kg)

Electrical Requirements (9000/9001)100-240V~, 50/60 Hz

9001 Only.....12 VDC @ 3 A max

Maximum Power Consumption65 watts

CPAP Pressure Range.....3-20 cm H₂O

Pressure Delay Time.....0, 10, 20, 30 or 45 minutes

Delay Pressure Range.....3 cm H₂O to prescription pressure

Operating Temperature Range41°F to 104°F (5°C to 40°C)

Operating Humidity Range5% to 95% RH non-condensing

Operating Atmospheric ConditionsSea level to 9,000 feet

Storage & Transportation Temperature Range.-4°F to 140°F (-20°C to +60°C)

Storage & Transportation Humidity Range.....5% to 95% RH non-condensing

Maximum Limited Pressure28 cmH₂O

Sound Level.....30 dBA (@ 10 cmH₂O, 1 meter from CPAP)

Filter specsStandard Filter > 3 micron particles

Optional Fine Particle Filter > 0.3 micron particles

Max flow rates	Model 9000	Model 9001	
	83 L/m	73 L/m	@ 7 cmH ₂ O (1/3 max pressure)
	134 L/m	118 L/m	@ 13 cmH ₂ O (2/3 max pressure)
	163 L/m	164 L/m	@ 20 cmH ₂ O (max pressure)

Pressure Stability

Between maintenance intervals Model 9000 = ± 2 cmH₂O
 (maintenance interval = 1 year) Model 9001 = ± 1 cmH₂O

Model 9000 Pressure swings (cmH₂O peak to peak at 500mL tidal volume)

Breaths per minute	10	15	20
@ 7 cmH ₂ O	1.0	1.0	1.0
@ 13 cmH ₂ O	1.0	1.0	1.0
@ 20 cmH ₂ O	1.3	1.3	1.4

Model 9001 Pressure swings (cmH₂O peak to peak at 500mL tidal volume)

Breaths per minute	10	15	20
@ 7 cmH ₂ O	1.0	1.0	1.1
@ 13 cmH ₂ O	1.0	1.0	1.1
@ 20 cmH ₂ O	1.1	1.1	1.2

Warranty.....Two-year limited

Equipment classification with respect to protection from electric shock.....Class II

Degree of protection from electric shock.....Type BF Applied Part

Degree of protection against ingress of liquids.....None

Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

Mode of operationContinuous

Definitions

Display – purpose, type, information displayed. The LCD display shows the following information in each model:

Model 9000	Model 9001	Description
Set Pres	Set Pres	CPAP pressure prescribed for the patient in cmH ₂ O
Altitude		Altitude setting one of four ranges
Fine Adj		Fine adjustment for fine tuning pressure using an external manometer
Dly Time	Dly Time	Delay time prescribed for the patient in minutes
Dly Pres	Dly Pres	Delay pressure prescribed for the patient in cmH ₂ O
Compl Hrs	Compl Hrs	Compliance meter, shows accumulated time “while breathing”
Unit Hrs	Unit Hrs	Unit hour meter, shows accumulated time of blower operation
Firmware	Firmware	Firmware version installed in the CPAP
Serial #	Serial #	Serial number of the CPAP
Language	Language	Language selection for LCD display messages
Off	Off	Blower is OFF AC or DC power is connected
Error	Error	An error has been detected by the firmware, the CPAP is in a fail safe mode

Pressure sensor – *purpose, type, sensing position, range, conditions of measurement*

The pressure sensor in Model 9001 measures the CPAP pressure inside the CPAP during operation. This measurement is used by the blower control circuit to correct the output pressure for mask leaks and altitude changes. The measured pressure is viewable by the PDA Remote. The pressure sensor is a silicon diaphragm type sensor with an input range of 0 to 25.4 cmH₂O. The measurement type of the sensor is “gauge” (relative to atmospheric pressure, ATPD).

General Info

Filters

Standard filter is designed to filter particles larger than 3 microns in size. The CPAP unit must have the standard filter in place during operation. This filter should be replaced if every 6 months or if it is damaged.

Optional Fine Particle Filter is designed to filter particles larger than 0.3 microns in size. The fine particle filter is placed in the filter area of the CPAP air intake with the standard filter placed over the fine particle filter. The optional fine particle filter should be replaced every 30 days or if it appears damaged.

Interdependence of controls – (settings that affect one another) The following settings on Model 9000 affect each other and must be set in the order listed to obtain the correct CPAP pressure:

- Set Pressure
- Altitude
- Fine Adjust

The CPAP Pressure Setting (Set Pres) determines the pressure for normal operation. The Set Pres value does not change when adjustments are made to Altitude or Fine Adjust settings.

The Manual Altitude Adjustment in Model 9000 allows a patient to adjust the pressure output of the CPAP to compensate for altitude changes when travelling. Changing the Altitude setting changes the CPAP pressure by up to 15% of the CPAP setting to correct for changes in barometric pressure at different altitudes.

The Fine Adjustment Control, available only during setup of the 9000 CPAP, adds or subtracts from the CPAP setting to obtain the desired pressure as measured by an external manometer.

The Delay Time and Delay Pressure operate the CPAP at a lower pressure for a preset period of time before reaching the prescribed CPAP pressure. The Altitude and Fine Adjust settings affect the Delay Pressure in the same way they affect the CPAP pressure.

Accuracy and Range of Displayed Values and Calibrated Controls Model 9001

Pressure sensor – max zero error = ± 1 cmH₂O

Sensitivity = $\pm 2\%$ FSO+4% reading

9000 Series CPAP Operation After Interruption and Restart of Power Supply

The 9000 Series CPAP will turn off if power is turned off or disconnected. When power is restored, the CPAP will return to the prescribed pressure setting within 2 seconds.

9000/9001 Series CPAP responds to the following conditions:

	Meets Spec:	9000	9001
Amb Temperature	+5°C to +40°C	✓	✓
Amb Humidity	10% to 95% RH	✓	✓
Atmos. Pressure	700 mbar to 1060 mbar	✓	✓
AC supply voltage range	-15% to +10% of declared nominal (100-240V)	✓	✓
DC supply voltage range	-15% to +25% of declared nominal (12 VDC)		✓

WARRANTY

DeVilbiss 9000/9001 Series nasal CPAP system is warranted to be free from defective workmanship and material for a period of two years from date of purchase. Any defective part(s) will be repaired or replaced at Sunrise Medical's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your Sunrise Medical provider or Sunrise Medical Service Department for instructions:

U.S.A. 800-333-4000, 814-443-4881

Canada 905-660-2459

Europe 44-138-444-6688.

NOTE—This warranty does not cover providing a loaner unit, compensating for costs incurred in rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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©Sunrise Medical 08.01
Sunrise Medical Ltd.
Sunrise Medical Pty. Limited