



**Impact** BioMedical CORP.

[www.impactbiomedical.com](http://www.impactbiomedical.com)

**1-800-327-1823**



## Respiratory

- Oxygen Concentrators
- Liquid Oxygen Systems
- H.E.P.A Filtration
- Foam Products
- Repair Parts
- Accessories

oxygen oxygen



*Authorized Warranty and  
Non Warranty Repair Center*



AUTHORIZED  
WARRANTY  
CENTER

WE PURCHASE  
OXYGEN EQUIPMENT

## RESPIRATORY EQUIPMENT SERVICE AND REPAIR DEPARTMENT

*Impact BioMedical Corp.* has been an AUTHORIZED warranty and non warranty repair center since 1989. We service Invacare, Mallinckrodt Puritan Bennett, Sunrise Medical / Devilbiss, Respironics, Nidek, Sequal and Chad therapeutics. This allows us to provide unmatched value and service you can only get from Impact BioMedical Corp. With procedures outlining return authorization, estimates and approval company assets will be managed PROFESSIONALLY AT THE HIGHEST STANDARD.

### RESPIRATORY EQUIPMENT & OXYGEN SYSTEMS SALES, RENTALS & REPAIRS

#### BRANDS SERVICED:

- MALLINCKRODT PURITAN BENNETT
- DEVILBISS SUNRISE HHG
- INVACARE CORPORATION
- RESPIRONICS
- RESMED
- NIDEK MEDICAL
- CHAD THERAPEUTICS
- AIRSEP
- SEQUAL
- EVO

#### REPAIR SERVICES FOR:

- CPAP & BILEVEL  
ALL MODELS
- OXYGEN CONCENTRATORS  
ALL MODELS
- LIQUID OXYGEN SYSTEMS  
BASES AND PORTABLES  
ALL MODELS
- SUCTION MACHINES  
ALL MODELS (BATTERY OPERATED ETC.)
- REGULATORS - OXYGEN  
ALL MODELS
- ANALYZERS - OXYGEN  
ALL MODELS
- 50 PSI COMPRESSORS  
ALL MODELS



1-800-327-1823

**OXYGEN  
CONCENTRATOR  
MONTHLY RENTAL  
AVAILABLE**

- NO REPAIR COST
- EXCHANGE WHEN NEEDED
- FREE FILTERS FOR RENTAL UNITS
- FIXED COST
- STOCK ONLY THE INVENTORY YOU NEED



**RECONDITIONED  
ECCENTRIC BEARING  
EXCHANGE ALL  
MODELS**



- COMPRESSORS AND ECCENTRIC BEARINGS
- ORDER BY COMPRESSOR TYPE AND CONCENTRATOR MODEL
- REQUIRES RETURN AUTHORIZATION

**RECONDITIONED  
COMPRESSOR  
EXCHANGE ALL  
MODELS**



Impact BioMedical, Corp.  
2408 Timberloch Place Suite B-3  
Phone (281) 292-6927 / FAX (281) 292-6911  
sales@impactbiomedical.com

TERMS: NET 30 ON OPEN ACCOUNT FROM DATE OF INVOICE  
WARRANTY: IMPACT BIOMEDICAL, INC.-MAKES NO WARRANTY OF ANY KIND, EXPRESSED OR IMPLIED, IN FACT OR BY LAW, EXCEPT THAT THE PRODUCTS SOLD IN THIS CATALOG WILL BE OF FINE QUALITY AND THE BUYER ASSUMES ALL RISK AND LIABILITY RESULTING FROM THE USE OF THESE PRODUCTS WHETHER USED SINGLY OR COMBINED WITH OTHER PRODUCTS OR SUBSTANCES.

**ALL PRICING SUBJECT TO CHANGE WITHOUT NOTICE!**

**1-800-327-1823**

**RETURN AUTHORIZATION REQUEST**

Impact BioMedical Corp.  
 2408 Timberloch Place  
 Suite B-3 The Woodlands, TX  
 77380

**ATTN:**

**FAX# OR E-MAIL**

**INSTRUCTIONS:** YOU HAVE REQUESTED SERVICE FOR YOUR EQUIPMENT.

WE NEED INFORMATION TO PROVIDE YOU A RETURN AUTHORIZATION NO.

NOTE: THE FOLLOWING ITEMS MAY NOT BE COVERED UNDER WARRANTY AND NON WARRANTY REPAIR SERVICES

- 1) SHIPPING EXPENSES
- 2) BOXES
- 3) CLEANING
- 4) FILTERS
- 5) ESTIMATE FEE

A) ALL REPAIR SERVICES WILL BE APPROVED BY YOUR OFFICE BEFORE ANY CHARGES APPLY.

B) PLEASE ISSUE A PURCHASE ORDER NUMBER IN THE SPACE PROVIDED TO ACKNOWLEDGE POSSIBLE CHARGES EXPLAINED ON 1 THRU 5 FOR WARRANTY REPAIRS AND NON-WARRANTY REPAIRS

C) PLEASE FILL OUT THE INFORMATION ON THIS FORM AND FAX BACK TO US AT (FAX#) 281-821-5635 (PHONE#) 281-821-5484

D) NOTE: THERE IS A \$50.00 ESTIMATE FEE ON ALL EQUIPMENT NOT APPROVED FOR SERVICE (PRICING ESTIMATES NEED INSTRUCTIONS RETURNED WITHIN 15 DAYS OR ESTIMATE FEE WILL BE INVOICED)

**NOTE!!!!**

- 1) ALL EQUIPMENT IS REQUIRED TO HAVE A RETURN AUTHORIZATION FROM IMPACT BIOMEDICAL CORP.
- 2) ALL EQUIPMENT RETURNED IN NEW BOXES
- 3) ALL EQUIPMENT FILTERS ARE REPLACED
- 4) ALL EQUIPMENT LEFT WITH OUT INSTRUCTIONS FOR 45 DAYS WILL BE DISPOSED OF.

YOUR COMPANY NAME:

PLEASE SHIP ME BOXES (CIRCLE YES OR NO)

PLEASE CHANGE FILTERS (CIRCLE YES OR NO)

PLEASE PICK UP EQUIP. (CIRCLE YES OR NO)

P.O. NO.

CONTACT FOR REPAIRS

CUST SIGNATURE (ACKNOWLEDGEMENT)

**PLEASE FAX BACK TO: 281-292-6927**

DESCRIPTION	HOURS IN	SERIAL NO.

# SAMPLE PRICING ESTIMATE

## Impact BioMedical Corp.

2408 Timberloch Place  
Suite B-3  
The Woodlands, TX 77380

CUST. NAME

ANY RESPIRATORY EQUIPMENT CO.  
12345 REPAIR STREET  
ANY PLACE, USA

### PLEASE FILL OUT COMPLETELY

INSTRUCTIONS: PLEASE FILL OUT 1 THRU 4 AND  
FAX TO 281-292-6911

THE FOLLOWING EQUIPMENT HAS BEEN SENT TO US FOR SERVICE AND  
REQUIRES INSTRUCTIONS. WE NEED THE FOLLOWING COMPLETED  
IMMEDIATELY AND FAXED TO 281-821-5635 TO PROCESS YOUR ORDER:

ATTN: YOUR CONTACT

FAX # 281-821-5484

DATE 06/25/2003

MODEL # PB590 WITH OCI

SERIAL # 0597G188708

WRTY STATUS SIEVE BEDS ONLY

HOURS IN 23809

SHIP TO LOCA- ANY PLACE

- 1) P.O. # \_\_\_\_\_
- 2) REPAIR : Y N (CIRCLE ONE) BILL CREDIT CARD Y OR N
- 3) HOW WOULD YOU LIKE THIS RETURNED? (CIRCLE ONE)
  - A) OUR COMPANY WILL PICK UP
  - B) SHIP UNIT
  - C) DISPOSE OF UNIT
- 4) SIGNATURE \_\_\_\_\_

### NOTICE!!!

ALL EQUIPMENT THAT DOES NOT RECEIVE INSTRUCTIONS ON SERVICE  
DESIRED WITH IN 15 BUSINESS DAYS WILL HAVE AN ESTIMSTE FEE INVOICED

ITEM	DESCRIPTION	QTY	WRTY	COST	TOTAL
DIAG	DIAGNOSTICS - WORK ORDER NO.12345 TECH TO PERFORM INSPECTION/TEST #1) FILTERS #2) COMPRESSOR #3) VIBRATION MOUNTS #4) TEST CYCLE PRESSURE #5) OUTLET PRESSURE #6) VALVING #7) LEAK	1	N	\$0.00	\$0.00
INCIDENTAL	INCIDENTAL	1	N	\$20.00	\$20.00
B-F010	INTAKE FILTER	1	N	\$3.75	\$3.75
C-F10	FOAM CABINET FILTER	1	N	\$0.50	\$0.50
F-F4	FELT PRE FILTER WITHOUT HOLE	2	N	\$0.32	\$0.64
F-F5	FELT PRE FILTER WITH HOLE	2	N	\$0.32	\$0.64
GR2V949	9 VOLT BATTERY	1	N	\$2.96	\$2.96
PW079	O-RING SIEVE BEDS	1	Y	\$0.00	\$0.00
PW095	18X18X28 BOX	1	Y	\$0.00	\$0.00
SHIP	SHIPPING AND HANDLING	1	N	\$20.00	\$20.00
WL	WARRANTY LABOR	1	Y	\$0.00	\$0.00

PLEASE FAX BACK TO 281-292-

TOTAL \$ 48.49

# Table of Contents

- 2) REPAIR INFORMATION
- 3) RENTAL CONCENTRATOR, COMPRESSOR REBUILD SERVICE, TERMS AND CONTACT INFORMATION
- 4) RETURN AUTHORIZATION FORM
- 5) SAMPLE PRICING ESTIMATE
- 6) TABLE OF CONTENTS
- 7) HEPA COMPRESSOR INTAKE FILTERS AND FINAL BACTERIA FILTER FOR CONCENTRATORS
- 8) HEPA COMPRESSOR INTAKE FILTERS AND PRE FILTER HOUSINGD
- 9) PRE FELT AND FOAM INTAKE FILTERS
- 10) PRE CUT FOAM CONCENTRATOR CABINET FILTERS
- 11) PRE CUT FOAM CONCENTRATOR CABINET FILTERS
- 12) HEPA COMPRESSOR INTAKE FILTER ASSYMBLIES
- 13) CONCENTRATOR FILTER SERVICE PACKS
- 14) CONCENTRATOR FILTER SERVICE PACKS
- 15) CONCENTRATOR FILTER SERVICE PACKS AND HEPA CAPSULE FILTERS
- 16) ASPIRATOR (SUCTION MACHINE) HYDROPHOBIC FILTERS AND NEB COMPRESSOR FILTERS
- 17) HUMIDIFIER ADAPTERS, FLOW KNOBS, TUBING
- 18) 2650 THOMAS STYLE COMPRESSOR PARTS
- 19) 100 THOMAS STYLE COMPRESSOR PARTS
- 20) GAST AND DUR-AIR STYLE COMPRESSOR SERVICE PARTS
- 21) 2619 THOMAS STYLE COMPRESSOR PARTS
- 22) 2618/17 THOMAS STYLE COMPRESSOR PARTS
- 23) MOTOR MOUNTS, DIAPHRAGMS, COOLING FANS AND CAPACITORS
- 24) REGULATOR SEALS, CYLINDER WRENCHES AND DISS CONNECTORS (CHRISTMAS TREE)
- 25) FAX ORDER FORM
- 26) CREDIT APPLICATION
- 27) EQUIPMENT RECONDITIONING POINT SHEET



**DESCRIPTION: HEPA INTAKE BACTERIA FILTER**  
**PART NO. B-FO1**

USE-OXYGEN CONCENTRATOR INLET FILTER  
 INLET CONNECTION-22 mm I.D.  
 OUTLET CONNECTION-22 mm O.D.

Price \$3.75 each  
 Pkg. Qty. 10

**FOR USE WITH OXYGEN CONCENTRATORS**



**DESCRIPTION: FINAL BACTERIA FILTER**  
**PART NO. B-FO2**

USE-OXYGEN CONCENTRATOR FINAL OUTLET BACTERIA FILTER  
 INLET CONNECTION-3/16"- 1/4" STEP BARB.  
 OUTLET CONNECTION- 3/16 - 1/4" STEP BARB

Price \$2.45 each  
 Pkg. Qty. 10



**DESCRIPTION: FINAL BACTERIA FILTER**  
**PART NO. B-FO3**

USE-OXYGEN CONCENTRATOR FINAL OUTLET BACTERIA FILTER  
 INLET CONNECTION-5/16" HOSE BARB.  
 OUTLET CONNECTION- 5/16 HOSE BARB

Price \$2.45 each  
 Pkg. Qty. 10

**FOR USE WITH OXYGEN CONCENTRATORS**



**DESCRIPTION: HEPA INTAKE BACTERIA FILTER**  
**PART NO. B-FO7**

USE-OXYGEN CONCENTRATOR COMPRESSOR INLET FILTER  
 INLET CONNECTION-22 mm I.D. / 3/8" FEMALE N.P.T. ADAPTER  
 OUTLET CONNECTION-22 mm O.D.



Pkg. Qty. 10

**1-800-327-1823**



**DESCRIPTION: EXTENDED LIFE INTAKE BACTERIA FILTER**

**PART NO. B-FO12**

*USE-OXYGEN CONCENTRATOR INLET FILTER  
INLET CONNECTION-22 mm*

**Price \$11.88 each**  
Pkg. Qty. 1



**DESCRIPTION: INLET FILTER**

**PART NO. B-FO9**

*USE-OXYGEN CONCENTRATOR INLET FILTER  
INLET CONNECTION-3/8" FEMALE N.P.T.*

**Price \$9.60 each**  
Pkg. Qty. 1

**FOR USE WITH OXYGEN CONCENTRATORS**



**DESCRIPTION: ROUND FELT PRE FILTER HOUSING**

**PART NO. P-FH1**

*USE -OXYGEN CONCENTRATOR COMPRESSOR INTAKE  
HEPA FILTER PRE FILTRATION*

**Price \$1.60 each**  
kg. Qty. 1



**DESCRIPTION: SQUARE FELT PRE FILTER HOUSING**

**PART NO. P-FH2**

*USE-OXYGEN CONCENTRATOR COMPRESSOR INTAKE  
HEPA FILTER PRE FILTRATION*

**Price \$2.75 each**  
Pkg. Qty. 1



**DESCRIPTION: FELT PRE FILTER WITH HOLE  
EXTENDED LIFE**

**PART NO. F-F5**

*USE -OXYGEN CONCENTRATOR PRE FILTRATION  
EXTENDS HEPA LIFE AND RESONATES INTAKE*

**Price .32¢ each**  
Pkg. Qty. 50



**DESCRIPTION: FELT PRE FILTER  
EXTENDED LIFE**

**PART NO. F-F4**

*USE -OXYGEN CONCENTRATOR PRE FILTRATION  
EXTENDS HEPA LIFE AND RESONATES INTAKE*

**Price .32¢ each**  
Pkg. Qty. 50



**DESCRIPTION: ROUND FELT PRE FILTER 3/8"**

**PART NO. F-F3**

*USE -OXYGEN CONCENTRATOR PRE FILTRATION  
EXTENDS HEPA LIFE AND RESONATES INTAKE*

**Price .32¢ each**  
Pkg. Qty. 50



**DESCRIPTION: FOAM PRE FILTER  
(RESONATOR)**

**PART NO. B-R1**

*USE -OXYGEN CONCENTRATOR PRE FILTRATION  
EXTENDS HEPA LIFE AND RESONATES INTAKE*

**Price .60¢ each**  
Pkg. Qty. 10



**DESCRIPTION: FOAM CABINET FILTER**  
**PART NO.C-F10**  
*PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE OF (PB) MALLINKCRODT CONCENTRATOR*

**Price .60¢ each Pkg. Qty. 10**



**DESCRIPTION: FOAM CABINET FILTER**  
**PART NO. C-F5**  
*PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE OF ALLIANCE CONCENTRATOR*

**Price .60¢ each Pkg. Qty. 10**



**DESCRIPTION: FOAM CABINET FILTER**  
**PART NO. C-F6**  
*PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE INVACARE CONCENTRATOR*

**Price .60¢ each Pkg. Qty. 10**



**DESCRIPTION: FOAM CABINET FILTER**  
**PART NO. C-F9**  
*PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE INVACARE PLATINUM CONCENTRATOR*

**Price .60¢ each Pkg. Qty. 10**



**DESCRIPTION: FOAM CABINET FILTER**  
**PART NO. C-F1**  
*PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE AIRSEP CONCENTRATOR*

**Price .60¢ each Pkg. Qty. 10**



**DESCRIPTION: FOAM CABINET FILTER**  
**PART NO. C-F14**  
*PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE MILLENNIUM CONCENTRATOR*

**Price .60¢ each Pkg. Qty. 10**



DESCRIPTION: FOAM CABINET FILTER

**PART NO. C-F11**

PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE  
*SOLARIS 5* CONCENTRATOR

Price .60¢ each Pkg. Qty. 10



DESCRIPTION: FOAM CABINET FILTER

**PART NO. C-F12**

PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE  
*SOLARIS 3* CONCENTRATOR

Price .60¢ each Pkg. Qty. 10



DESCRIPTION: FOAM CABINET  
FILTER

**PART NO. C-F3**

PRE CUT FILTER FOR OUTSIDE CABINET  
ON INTAKE DEVILBISS *MC44, MC64,*  
*MC84, PVO2,* CONCENTRATOR

Price .60¢ each Pkg. Qty. 10



DESCRIPTION: FOAM CABINET FILTER

**PART NO. C-F18**

PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE  
*SEQUAL* CONCENTRATOR

Price .60¢ each Pkg. Qty. 10



DESCRIPTION: FOAM CABINET FILTER

**PART NO. C-F4**

PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE OF  
*NIDEK NUVO* CONCENTRATOR

Price .60¢ each Pkg. Qty. 10



DESCRIPTION: FOAM CABINET FILTER

**PART NO. C-F7**


PRE CUT FILTER FOR OUTSIDE CABINET ON INTAKE  
*NIDEK MARK 5* CONCENTRATOR

Price \$5.40 each Pkg. Qty. 10



**DESCRIPTION: 6 PIECE ASSEMBLY- INTAKE BACTERIA FILTER W/EXTENDED LIFE PRE FILTERS AND HOUSING**  
**PART NO. i-A1**

Price \$ 7.78 each Pkg. Qty. 1




**DESCRIPTION: 3 PIECE ASSEMBLY-INTAKE BACTERIA FILTER W/ROUND PRE FILTER AND HOUSING**  
**PART NO. i-A2**

Price \$ 5.67 each Pkg. Qty. 1



**DESCRIPTION: 2 PIECE ASSEMBLY-INTAKE BACTERIA FILTER W/FOAM PRE FILTER**  
**PART NO. i-A3**

Price \$ 4.35 each Pkg. Qty. 1



**DESCRIPTION: 5 PIECE ASSEMBLY-EXTENDED LIFE PRE FILTERS AND HOUSING**  
**PART NO. i-PFA1**

Price \$ 4.03 each Pkg. Qty. 1



DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P7**

*SOLARIS 5 CONCENTRATOR*

Price \$ 4.67 each  
Pkg. Qty. 1



DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P6**

*SOLARIS 5 CONCENTRATOR*

Price \$ 4.67 each  
Pkg. Qty. 1



DESCRIPTION: 3 PIECE FILTER SERVICE PACK (

**PART NO. S-P9**

*SOLARIS 3 CONCENTRATOR*

Price \$ 4.67 each  
Pkg. Qty. 1



DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P8**

*SOLARIS 3 CONCENTRATOR*

Price \$ 4.67 each  
Pkg. Qty. 1



DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P11**

*DEVILBISS PVO2-MC44-MC64-MC84*

Price \$ 4.67 each  
Pkg. Qty. 1




DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P12**

*DEVILBISS PVO2-MC44-MC64-MC84*

Price \$ 4.67 each  
Pkg. Qty. 1






DESCRIPTION: 6 PIECE FILTER SERVICE PACK

**PART NO. S-P1**

*PB CONCENTRATOR*



Price \$ 5.63 each  
Pkg. Qty. 1




DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P2**

*PB CONCENTRATOR*



Price \$ 4.67 each  
Pkg. Qty. 1



DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P3**

*ALLIANCE CONCENTRATOR*



Price \$ 4.95 each  
Pkg. Qty. 1




DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P13**

*MILLENNIUM CONCENTRATOR*



Price \$ 4.95 each  
Pkg. Qty. 1




DESCRIPTION: 3 PIECE FILTER SERVICE PACK

**PART NO. S-P5**

*H300 & BX5000 CONCENTRATOR*



Price \$ 12.48 each




DESCRIPTION: 4 PIECE FILTER SERVICE PACK

**PART NO. S-P4**

*INVACARE CONCENTRATOR*




Price \$ 8.35 each  
Pkg. Qty. 1



**DESCRIPTION: 4 PIECE FILTER SERVICE PACK**  
**PART NO. S-P10**  
 INVACARE *IRC5LX* CONCENTRATOR



**Price \$ 13.40 each**  
 Pkg. Qty. 1



**DESCRIPTION: 4 PIECE FILTER SERVICE PACK**  
**PART NO. S-P11**  
 INVACARE *IRC5LX PLATINUM* CONCENTRATOR




**Price \$ 13.40 each**  
 Pkg. Qty. 1



**DESCRIPTION: 1 PIECE ELEMENT FILTER**  
**PART NO. NT-002**  
 NIDEK CONCENTRATOR

**Price \$ 6.58 each**  
 Pkg. Qty. 1



**DESCRIPTION: 2 PIECE FILTER SERVICE PACK**  
**PART NO. S-P16**  
 NIDEK CONCENTRATOR



**Price \$ 13.06 each**  
 Pkg. Qty. 1



**Price \$ 16.50 each**  
 Pkg. Qty. 1

**DESCRIPTION: 6" HEPA CAPSULE FILTER**  
**PART NO. i-H3**  
 USE - OXYGEN CONCENTRATOR INTERNAL IN LINE FILTRATION  
 INLET CONNECTION - 3/8" FEMALE THREADED  
 OUTLET CONNECTION - 3/8" FEMALE THREADED  
 EFFICIENCY - >99.97 @ 0.3 MICRON  
 MAXIMUM OPERATING TEMPERATURE - 212F @ 0 PSI



**Price \$ 16.50 each**  
 Pkg. Qty. 1

**DESCRIPTION: 3" HEPA CAPSULE FILTER**  
**PART NO. i-H2**  
 USE - OXYGEN CONCENTRATOR INTERNAL IN LINE FILTRATION  
 INLET CONNECTION - 3/8" FEMALE THREADED  
 OUTLET CONNECTION - 3/8" FEMALE THREADED  
 EFFICIENCY - >99.97 @ 0.3 MICRON  
 MAXIMUM OPERATING TEMPERATURE - 212F @ 0 PSI



**DESCRIPTION: HYDROPHOBIC BACTERIA FILTER**  
**PART NO. B-FO4**

USE - ASPIRATORS (SUCTION MACHINES) VACUUM INLET

INLET CONNECTION - 3/8" THREADED MALE NPT  
OUTLET CONNECTION - 3/16" - 1/4" STEP BARB

Price \$ 2.60 each  
Pkg. Qty. 10



**DESCRIPTION: HYDROPHOBIC BACTERIA FILTER**  
**PART NO. B-FO5**

USE - ASPIRATORS (SUCTION MACHINES) VACUUM INLET

INLET CONNECTION - 3/8" THREADED MALE NPT  
OUTLET CONNECTION - 3/4" - 5/16" STEP BARB

Price \$ 2.60 each  
Pkg. Qty. 10



**DESCRIPTION: HYDROPHOBIC BACTERIA FILTER**  
**PART NO. B-FO6**

USE - ASPIRATORS (SUCTION MACHINES) VACUUM INLET

INLET CONNECTION - 3/16" STEP BARB  
OUTLET CONNECTION - 3/16" - STEP BARB

Price \$ 2.60 each  
Pkg. Qty. 10



**DESCRIPTION: HYDROPHOBIC BACTERIA FILTER**  
**PART NO. B-FO8**

USE - ASPIRATORS (SUCTION MACHINES) VACUUM INLET

INLET CONNECTION - 3/8" THREADED MALE NPT  
OUTLET CONNECTION - 3/16" - STEP BARB

Price \$ 2.60 each  
Pkg. Qty. 10

DESCRIPTION: NEBULIZER COMPRESSOR INTAKE  
FILTER INSPIRATION 929 & 626

**PART NO. N-929F**

Price \$ .58¢ each

Pkg. Qty. 50



DESCRIPTION: NEBULIZER COMPRESSOR INTAKE  
FILTER FOR INSPIRATION 323

**PART NO. N-323F**

Price \$ 1.14 each

Pkg. Qty. 50



DESCRIPTION: NEBULIZER COMPRESSOR INTAKE  
FILTER FOR PULMOMATE

**PART NO. N-PM1**

Price \$ .35¢ each

Pkg. Qty. 50



DESCRIPTION: NEBULIZER COMPRESSOR INTAKE  
FILTER FOR PULMOAIDE

**PART NO. N-PM2**

Price \$ .5¢ each

Pkg. Qty. 50



DESCRIPTION: NEBULIZER COMPRESSOR INTAKE  
FILTER FOR AEROMAX, SPORT NEB, INVACARE  
PASSPORT & SCHUCO

**PART NO. N-AF1**

Price \$ .10¢ each

Pkg. Qty. 50





**DESCRIPTION: CLEAR 6" TUBE WITH PLASTIC HUMIDIFIER ADAPTER**  
**PART NO. H-A6**  
 USE - OXYGEN CONCENTRATOR HUMIDIFIER ADAPTER

**Price \$ 2.50 each Pkg. Qty. 10**



**DESCRIPTION: 2" CHROME HUMIDIFIER ADAPTER**  
**PART NO. H-A2**  
 USE - HUMIDIFIER ADAPTER

**Price \$4.50 each Pkg. Qty. 10**



**DESCRIPTION: CLEAR 3" TUBE WITH PLASTIC HUMIDIFIER ADAPTER**  
**PART NO. H-A4**  
 USE - OXYGEN CONCENTRATOR HUMIDIFIER ADAPTER

**Price \$ 2.50 each Pkg. Qty. 10**



**DESCRIPTION: PLASTIC HUMIDIFIER ADAPTER**  
**PART NO. H-PB21**  
 USE - OXYGEN CONCENTRATOR HUMIDIFIER ADAPTER

**Price \$ 4.50 each Pkg. Qty. 10**



**DESCRIPTION: FLOW METER KNOB**  
**PART NO. F-MK1**

**Price \$ 3.20 each Pkg. Qty. 10**



**DESCRIPTION: PLASTIC HUMIDIFIER ADAPTER WITH OUT HOSE**  
**PART NO. IT025**

**Price \$ 1.50 each Pkg. Qty. 10**



**DESCRIPTION: OXYGEN OUTLET SPOUT**  
**PART NO. H-PB4**  
 USE - OXYGEN CONCENTRATOR OUTLET FOR HUMIDIFIER ADAPTER

**Price \$ 4.80 each Pkg. Qty. 1**



**DESCRIPTION: OXYGEN OUTLET INVACARE**  
**PART NO. H-INV8**  
 USE - OXYGEN CONCENTRATOR OUTLET FOR HUMIDIFIER ADAPTER

**Price \$ 3.75 each Pkg. Qty. 10**

**DESCRIPTION: VINYL TUBING 1/2" O.D. X 1/4" I.D. 1/8" WALL X 100' ROLL**  
**PART NO. RHR-001**  
 USE - OXYGEN CONCENTRATOR OUTLET FOR HUMIDIFIER ADAPTER CONNECTIONS

**Price \$ 45.00 each Pkg. Qty. 100 FT. ROLL**



# 2650 THOMAS STYLE COMPRESSOR SERVICE PARTS



**DESCRIPTION: 2650 STYLE COMPRESSOR TOP END REBUILD KIT**

## PART NO. T-2650

USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
KIT INCLUDES-2 PISTON CUPS, 2 "H" GASKATS, 4 FLAPPER VALVES, 2 SLEEVE O-RINGS, 2 RESTRAINTS

Price \$ 25.20 each  
Pkg. Qty. 1 KIT



**DESCRIPTION: 2650 STYLE COMPRESSOR PISTON CUP SEAL**

## PART NO. T-18

USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
SOLD AS EACH

Price \$ 5.05 each  
Pkg. Qty. 1



**DESCRIPTION: 2650 STYLE COMPRESSOR FLAPPER VALVE**

### PART NO. T-2

Price \$ 1.40 each



**DESCRIPTION: 2650 STYLE COMPRESSOR "H" GASKET**

### PART NO. T-19

Price \$ 2.96 each



**DESCRIPTION: 2650 STYLE COMPRESSOR FLAPPER VALVE RESTRAINT**

### PART NO. T-6

Price .55¢ each



**DESCRIPTION: 2650 STYLE COMPRESSOR PISTON SLEEVE O-RING**

### PART NO. T-16

Price \$ 2 .60 each



**DESCRIPTION: 2650 STYLE COMPRESSOR PISTON SLEEVE**

### PART NO.

Price \$ 7.00 each



**DESCRIPTION: 2650 STYLE COMPRESSOR ECCENTRIC BEARING**

### PART NO. T-8

Price \$ 5.60 each



# *100 THOMAS STYLE COMPRESSOR SERVICE PARTS*



**DESCRIPTION: 100 STYLE COMPRESSOR TOP END REBUILD KIT**

## **PART NO. T-100**

*USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
KIT INCLUDES- 1 PISTON CUP, 1 HEAD GASKET, 3 FLAPPER VALVES, 1 CYLINDER SLEEVE, 1 O-RING FILTER COVER, 1 O-RING CYLINDER SLEEVE, 1 FELT FILTER, 1 GASKET RETAINER, 2 VALVE STOP, 2 FLAPPER VALVE SCREW, 4 RETAINER SCREW, 1 RETAINER ASSEMBLY, 1 VALVE KEEPER, 1 IN LINE MUFFLER*



**Price \$ 36.75 each**  
**Pkg. Qty. 1**

<p><b>DESCRIPTION:</b> <b>100 STYLE COMPRESSOR PISTON SLEEVE O-RING</b></p>  <p><b>PART NO. T-25</b></p> <p><b>Price \$ 2.90 each Pkg. Qty. 1</b></p>	<p><b>DESCRIPTION: 100 STYLE COMPRESSOR PISTON CUP SEAL</b></p>  <p><b>PART NO. T-29</b></p> <p><b>Price \$ 9.30 each Pkg. Qty. 1</b></p>
<p><b>DESCRIPTION:</b> <b>100 STYLE COMPRESSOR HEAD O-RING</b></p>  <p><b>PART NO. T-26</b></p> <p><b>Price \$ 3.36 each Pkg. Qty. 1</b></p>	<p><b>DESCRIPTION: 100 STYLE COMPRESSOR PISTON SLEEVE</b></p>  <p><b>PART NO. T-30</b></p> <p><b>Price \$ 10.00 each</b></p>
<p><b>DESCRIPTION: 100 STYLE COMPRESSOR FELT INTAKE FILTER</b></p>  <p><b>PART NO. T-28</b></p> <p><b>Price .50 each Pkg. Qty. 1</b></p>	<p><b>DESCRIPTION: 100 STYLE COMPRESSOR INLINE MUFFLER</b></p>  <p><b>PART NO. T-27</b></p> <p><b>Price \$ 27.50 each Pkg. Qty. 1</b></p>

## GAST STYLE COMPRESSOR SERVICE PARTS



DESCRIPTION:  
GAST STYLE COMPRESSOR "H"  
GASKET



**PART NO. SWC002**

Price \$3.26 each

DESCRIPTION:  
GAST STYLE COMPRESSOR  
PISTON CUP SEAL



**PART NO. T-17**

Price \$ 9.50 each

DESCRIPTION: GAST STYLE  
COMPRESSOR PISTON SLEEVE O-  
RING



**PART NO. SWC001**

Price \$ 2.60 each

DESCRIPTION: GAST STYLE  
COMPRESSOR PISTON SLEEVE



**PART NO. T-32**

Price \$ 10.00

## DURA-AIR STYLE COMPRESSOR SERVICE PARTS



DESCRIPTION: DURA AIR STYLE COMPRESSOR TOP END REBUILD KIT  
**PART NO. MC-001**

USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
KIT INCLUDES-2 PISTON CUPS, 2 "H" GASKETS, 4 FLAPPER VALVES, 4 TUBE  
O-RINGS, 2 SLEEVE O-RINGS



Price \$ 33.33 each  
Pkg. Qty. 1

DESCRIPTION: DURA-AIR STYLE COMPRESSOR PISTON CUP SEAL  
**PART NO. T-31**

USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
SOLD AS EACH



Price \$ 7.70 each Pkg. Qty. 1

# 2619 STYLE THOMAS COMPRESSOR SERVICE PARTS



**DESCRIPTION: 2619 STYLE COMPRESSOR TOP END REBUILD KIT**

## PART NO. T-2619

USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
 KIT INCLUDES-2 PISTON CUPS, 2 "H" GASKATS, 4 FLAPPER VALVES, 4 TUBE ORINGS,  
 2 SLEEVE O-RINGS, 2 RESTRAINTS

Price \$ 15.00 each  
 Pkg. Qty. 1 KIT



**DESCRIPTION: 2619 STYLE COMPRESSOR PISTON CUP SEAL**

## PART NO. T-1

USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
 SOLD AS EACH

Price \$ 2.05 each  
 Pkg. Qty. 1



**DESCRIPTION: 2619 STYLE COMPRESSOR FLAPPER VALVE**

**PART NO. T-2** Price \$ 1.40 each



**DESCRIPTION: 2619 STYLE COMPRESSOR FLAPPER VALVE RESTRAINT**

**PART NO. T-6** Price .55¢ each



**DESCRIPTION: 2619 STYLE COMPRESSOR CROSSOVER TUBE O-RING**

**PART NO. T-5** Price .05¢ each



**DESCRIPTION: 2619 STYLE COMPRESSOR "H" GASKET**

**PART NO. T-3** Price \$ 1.50 each



**DESCRIPTION: 2619 STYLE COMPRESSOR PISTON SLEEVE O-RING**

**PART NO. T-4** Price .70¢ each



**DESCRIPTION: 2619 STYLE PISTON SLEEVE**

**PART NO. T-24**

Price \$ 7.00 each



**DESCRIPTION: 2619 STYLE COMPRESSOR ECCENTRIC BEARING**

**PART NO. T-8**

Price \$ 5.60 each



# 2617/18 THOMAS STYLE COMPRESSOR SERVICE PARTS



DESCRIPTION: 2617/18 STYLE COMPRESSOR TOP END REBUILD KIT

## PART NO. T-2617

USE - OXYGEN CONCENTRATOR COMPRESSOR MAINTENANCE:  
KIT INCLUDES-2 PISTON CUPS, 2 "H" GASKATS, 4 FLAPPER VALVES, 4 TUBE ORINGS,  
2 SLEEVE O-RINGS, 2 RESTRAINTS

Price \$ 23.08 each

Pkg. Qty. 1 KIT



DESCRIPTION:  
2617/18 STYLE COMPRESSOR PISTON CUP SEAL

## PART NO. T-12

Price \$ 4.05 each

Pkg. Qty. 1



DESCRIPTION: 2617/18 STYLE  
COMPRESSOR FLAPPER VALVE

PART NO. T-2 Price \$ 1.40 each



DESCRIPTION:  
2617/18 STYLE COMPRESSOR "H"  
GASKET

PART NO. T-13 Price \$ 2.70 each



DESCRIPTION: 2617/18 STYLE  
COMPRESSOR FLAPPER VALVE  
RESTRAINT

PART NO. T-6 Price .55¢ each



DESCRIPTION:  
2617/18 STYLE COMPRESSOR  
PISTON SLEEVE O-RING

PART NO. T-14 Price \$ 1.24 each



DESCRIPTION: 2617/18 STYLE  
COMPRESSOR CROSSOVER TUBE O-RING

PART NO. T-5 Price .05¢ each



DESCRIPTION:  
2617/18 STYLE PISTON SLEEVE

PART NO. T-34

Price \$ 10.00 each



DESCRIPTION:  
2617/18 STYLE COMPRESSOR  
ECCENTRIC BEARING

PART NO. T-8

Price \$ 5.60 each



DESCRIPTION: DIAMOND STYLE  
MOTOR MOUNT

**PART NO. M-M2**

USE-ISOLATES VIBRATION OF  
COMPRESSOR



Price \$ 2.50 each  
Pkg. Qty. 1

DESCRIPTION: SQUARE STYLE  
MOTOR MOUNT

**PART NO. M-M3**

USE-ISOLATES VIBRATION OF  
COMPRESSOR



Price \$ 3.50 each  
Pkg. Qty. 1

DESCRIPTION: POST STYLE MOTOR  
MOUNT

**PART NO. M-M1**

USE-ISOLATES VIBRATION OF  
COMPRESSOR



Price \$ 2.50 each  
Pkg. Qty. 1

DESCRIPTION: BUBBLE STYLE MOTOR  
MOUNT 6 LB. & 9 LB.

**PART NO. M-M4**

USE-ISOLATES VIBRATION OF  
COMPRESSOR



Price \$ 4.80 each  
Pkg. Qty. 1

DESCRIPTION: DIAMOND STYLE  
MOTOR MOUNT  
(SMALL)

**PART NO. M-M6**

USE-ISOLATES VIBRATION OF  
COMPRESSOR



Price \$ 3.50 each  
Pkg. Qty. 1

DESCRIPTION: FRONT SOLID  
DIAPHRAGM

**PART NO. GDP-1**

PB CONCENTRATOR



DESCRIPTION: FRONT SOLID  
DIAPHRAGM

**PART NO. GDA-1**

ALLIANCE CONCENTRATOR

Price \$ 3.50 each  
Pkg. Qty. 1

DESCRIPTION: REAR  
PERFERATED DIAPHRAGM

**PART NO. GDP-2**

PB CONCENTRATOR



Price \$ 3.50 each  
Pkg. Qty. 1

DESCRIPTION: REAR  
PERFERATED DIAPHRAGM

**PART NO. GDA-2**

ALLIANCE CONCENTRATOR

DESCRIPTION: COOLING FAN

**PART NO. T-FN**

USE-OXYGEN CONCENTRATOR  
COOLING FAN



Price \$23.00 each  
Pkg. Qty. 1

DESCRIPTION:  
CAPACITOR 7.5MF

**PART NO. R-C3**



Price \$ 7.85 each  
Pkg. Qty. 1

DESCRIPTION:  
CAPACITOR 15MF

**PART NO. R-C2**



Price \$19.50 each  
Pkg. Qty. 1

DESCRIPTION:CAPACIT  
OR 7MF

**PART NO. R-C1**



Price \$ 9.95 each  
PKG. QTY. 1



DESCRIPTION: TAMPER EVIDENT WRAP & LOCK DUST COVER

**PART NO. W-001**

Price .10¢ each  
Pkg. Qty. 50

COLORS:  
WHITE



DESCRIPTION: YOKE SEAL -PLASTIC

**PART NO. W-006**

COLORS:  
YELLOW  
GREEN

Price .10¢ each  
Pkg. Qty. 50



COLORS:  
GREEN  
RED  
BLACK

DESCRIPTION: ECONOMY CYLINDER POST VALVE WRENCH

**PART NO. W-003**

Price .45¢ each  
Pkg. Qty. 50



DESCRIPTION: YOKE SEAL -BRASS WITH RUBBER O-RING

**PART NO. W-RS2**

USE -OXYGEN REGULATOR SEAL

Price .45¢ each  
Pkg. Qty. 50



DESCRIPTION: ECONOMY CYLINDER POST VALVE WRENCH WITH CHAIN

**PART NO. W-007**

Price .60¢ each  
Pkg. Qty. 100

COLORS:  
GREEN  
RED  
BLACK



DESCRIPTION: GREEN PLASTIC DISS CONNECTOR

**PART NO. H-DC1**

USE- OXYGEN CONNECTOR-ETC.

Price .45¢ each  
Pkg. Qty. 50



**CREDIT APPLICATION**

COMPANY FULL LEGAL NAME: \_\_\_\_\_

ADDRESS: \_\_\_\_\_

PHONE# \_\_\_\_\_ FAX# \_\_\_\_\_

MAILING ADDRESS: \_\_\_\_\_

CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

YEARS IN BUSINESS \_\_\_\_\_ CREDIT REQUESTED \_\_\_\_\_ (U.S. DOLLARS)

BUSINESS STRUCTURE—CORPORATION ( ) / PARTNERSHIP ( ) / PROPRIETORSHIP ( ) / OTHER ( )

**COMPANY PRINCIPALS RESPONSIBLE FOR BUSINESS TRANSACTIONS**

NAME / TITLE

ADDRESS / CITY / STATE / ZIP


**BANK REFERENCES**

**INSTITUTION NAME      ACCOUNT NO.      CONTACT NAME      PHONE #      FAX #**


**TRADE REFERNCES**

**COMPANY NAME      CONTACT NAME      PHONE #      FAX #      ACCOUNT #**


**CONFIRMATION OF INFORMATION ACCURACY AND RELEASE OF AUTHORITY TO VERIFY**

I HEREBY CERTIFY THAT THE INFORMATION IN THIS ACCOUNT APPLICATION IS CORRECT. THE INFORMATION IN THIS ACCOUNT APPLICATION IS FOR USE BY IMPACT BIO-MEDICAL, INC. IN DETERMINING THE AMOUNT AND CONDITIONS OF CREDIT TO BE EXTENDED. I UNDERSTAND THAT TMR COMPANIES INC. MAY ALSO UTILIZE THE OTHER SOURCES OF CREDIT WHICH IT CONSIDERS NECESSARY IN MAKING THIS DETERMINATION. FURTHER, I HEREBY AUTHORIZE THE BANK AND TRADE REFERENCES LISTED IN THIS ACCOUNT APPLICATION TO RELEASE THE INFORMATION NECESSARY TO ASSIST IMPACT BIOMEDICAL, NC. IN ESTABLISHING A LINE OF CREDIT.

**SIGNATURE** \_\_\_\_\_ **DATE** \_\_\_\_\_ **TITTLE** \_\_\_\_\_

**PHONE # 281-292-6911  
FAX # 281-292-6927**

WE SELL ALL  
MODELS OF  
RECONDITIONED  
OXYGEN  
CONCENTRATORS,  
REGULATORS AND  
LIQUID OXYGEN  
SYSTEMS

!!!!!!  
**SAVE**  
!!!!!!  
RECONDITIONED  
EQUIPMENT FOR  
**SALE**

\$\$\$\$\$  
WE PURCHASE  
\$\$\$\$\$  
OXYGEN  
CONCENTRATORS,  
REGULATORS AND  
LIQUID OXYGEN  
SYSTEMS



### 1—QUALITY

To provide a safe and reliable product  
Impact BioMedical Corp. **does not “modify” (cutting  
top and refilling) sieve beds eliminating possible  
liability due to altering original factory design.**

### 2—LICENSED

Impact BioMedical Corp. is properly licensed to  
recondition respiratory equipment through the  
**TDH Bureau of Food and Drug Safety**

### 3—RECONDITIONING EXPERIENCE

We have serviced over 100,000 oxygen systems since  
1989. We recondition equipment with procedures and  
standards for the professional that requires quality and  
value.

### 4—WARRANTY

WARRANTY: 1 YEAR STANDARD-ADDITIONAL  
WARRANTY CAN BE PURCHASED/\*\*CALL FOR  
MODELS AND AVAILABILITY\*\*

CALL FOR  
AVAILABILITY &  
PRICING

1-800-327-1823

Impact Biomedical Corp.  
2408 Timberloch Place Suite B-3  
The Woodlands ,TX 77380  
Phone 800-327-1823 Fax 281-292-6927



## **WE STOCK A FULL LINE OF FILTERS AND REPLACEMENT PARTS**

- **HEPA INTAKE BACTERIA FILTERS**
- **FINAL BACTERIA FILTERS**
- **PRE FILTER HOUSINGS**
- **FOAM PRE FILTERS**
- **PRE CUT CABINET FILTERS**
- **INTAKE BACTERIA FILTER ASSEMBLIES**
- **CONCENTRATOR FILTER SERVICE PACKS**
- **INLINE HEPA FILTER CAPSULES**
- **SUCTION MACHINE (ASPIRATORS) HYDROPHOBIC BACTERIA FILTERS**
- **NEBULIZER COMPRESSOR FILTERS**
- **COMPRESSOR KITS AND PARTS**
- **COOLING FANS**
- **COMPRESSOR VIBRATION MOUNTS (MOTOR MOUNTS)**
- **MOTOR (COMPRESSOR) RUN CAPACITORS**
- **OXYGEN CONCENTRATOR HUMIDIFIER ADAPTERS AND CONNECTORS**
- **OXYGEN OUTLETS AND CYLINDER WRENCHES**
- **CYLINDER SEALS AND CAPS**
- **RECONDITIONED EQUIPMENT**
- **COMPRESSOR AND EQUIPMENT REPAIR SERVICES**
- **AUTHORIZED WARRANTY REPAIR**

**CUSTOM BOXES FOR ALL MODELS OF CONCENTRATORS \$25.88 EACH**